

## Cystatin C free serum

## **Product Summary**

Source:

Pooled normal human serum. Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

Applications:

Cystatin C free serum can be used as a matrix for standard and calibrator preparation.

Analysis:

According to ELISA testing cystatin C free serum does not contain more than 2 ng/ml of cystatin C.

Purification

Immunoaffinity chromatography of pooled normal human serum.

Gel for affinity column was prepared by coupling anti-cystatin C monoclonal antibodies with BrCNactivated Sepharose CL 4B.

Presentation:

Frozen liquid.

Storage Conditions:

-20 °C (between -15 and -30 °C allowed)

Material Safety Note:

This product is sold for research use only. Standard Laboratory Practices should be followed when handling this material.

Additional product information:

Cystatin C Assay Notes

References:

Search Cystatin C references from PubMed. See all Cystatin C references.

See other  $\underline{\text{cardiac markers}}$  and  $\underline{\text{blood coagulation}}$  products.

Related products:

Cat# 2-CYSC: Cystatin C MoAb
Cat# 6-CYS: Polyclonal anti-cystatin C
Cat# 8-CYC: Human cystatin C

Cat# 8-CYC-rh: Cystatin C, recombinant human

Price: \$1,480.00

Catalog Id: 11-CYC-fs

**Product Categories**: Sera & Plasma

Product Page: https://www.advimmuno.com/product/cystatin-c-free-serum/

**Product Attributes** 

- Sizes:: 50 ml

Advanced ImmunoChemical, Incorporated | advimmuno.com

562-434-4676 | tech@advimmuno.com | order@advimmuno.com 111 West Ocean Blvd, Landmark Square, 4th Floor, Long Beach, CA 90802