

Myeloperoxidase free serum

Product Summary

Source:

Pooled normal human serum.

Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

Applications:

Myeloperoxidase free serum can be used for preparation of human myeloperoxidase standards and calibrators.

Analysis:

No detectable amounts of human myeloperoxidase.

Purification:

Affinity chromatography of pooled normal human serum.

Gel for affinity column was prepared by coupling anti-myeloperoxidase monoclonal antibodies (several monoclonal antibodies with different epitope specificity) with BrCN-activated Sepharose CL 4B.

Presentation:

Frozen liquid.

Storage Conditions:

-20 °C (between -15 and -30 °C allowed)

Material Safety Note:

This product is sold for research use only. Standard Laboratory Practices should be followed when handling this material.

Additional product information:

[Myeloperoxidase \(MPO\) Assay Notes](#)

References:

[Search MPO references from PubMed. See all MPO references.](#)

Price: \$1,480.00

Catalog Id: 11-MPO-fs

Product Categories: [Sera & Plasma](#)

Product Page: <https://www.advimmuno.com/product/myeloperoxidase-free-serum/>

Product Attributes

- Sizes:: 50 ml

Advanced ImmunoChemical, Incorporated | [advimmuno.com](https://www.advimmuno.com)
562-434-4676 | tech@advimmuno.com | order@advimmuno.com
111 West Ocean Blvd, Landmark Square, 4th Floor, Long Beach, CA 90802